#### DRYING OVEN VERIFICATION PROCEDURE #2

Revised 10/12/1995

Equipment Checked: DRYING OVEN

#### Purpose:

This method provides instructions for verifying the dial settings on general purpose drying ovens.

### **Inspection Equipment Required:**

- 1. A calibrated thermometer graduated in 1.0 degress C increments having a range which includes the temperature range to be checked.
- 2. A brass thermometer well to retain heat while the oven door is open. This is essential for a constant temperature reading
- 3. A clothes pin to hold the thermometer in such a manner as to enable the operator to read the scale easily from outside or inside the oven.

#### Tolerance:

Drying ovens shall be capable of maintaining a constant temperature range listed in the appropriate test methods.

#### **Procedure:**

- 1. Place the thermometer inside the brad well with the clothes pion attached to the thermometer. Position the thermometer on the shelf where the samples are normally dried.
- 2. Take the first reading at least 1 hour after closing the oven (oven should remain undisturbed).
- 3. Take as many readings as necessary to determine if the temperature range is within the specified tolerance (three consecutive readings, taken no less then ½ hour apart within the tolerance allowed are adequate).
- 4. Adjust the temperature of the oven if an observed temperature reading is outside the tolerance specified (aloe at least ½ hour for the temperature to stabilize between each adjustment). Return to Step 3.

# **EQUIPMENT CALIBRATION RECORD**

## **Calibration Procedure No.2 Ovens**

Previous Calibration Date:			Calibrated by:
			Next Due:
Frequency: Indeper	ndent Assurance – 12 m	onths, Acceptance Testing	g – Prior to use on project
Action Recommend	ed:		
Repair	Replace	None Other	
Calibration Equipment:		Serial Number:	
Calibrated Thermon	neter		
MODEL		MODEL	
S/N		S/N	
SET @	ACTUAL	SET @	ACTUAL
MODEL		MODEL	
S/N		S/N	
SET @	ACTUAL	SET @	ACTUAL
MODEL		MODEL	
S/N		S/N	
SET @	ACTUAL	SET @	ACTUAL
MODEL		MODEL	
S/N		S/N	
SET @	ACTUAL	SET @	ACTUAL
MODEL		MODEL	
MODEL		MODEL	
S/N	4.000	S/N	A 07115:
SET @	ACTUAL	SET @	ACTUAL